

# Help

**Data format:** folder containing only recordings of the same machine and same protocol (can lead to text display “format error” or “shape error”)

To load faster (from ~20s to 1-2s), copy the folder ‘foxexe’ from T/studies on your computer, and go to ‘forexe/dist/app/app.exe’ to create a new shortcut

## How to use the interface:

Click on select (top left menu), select the machine, a window open to select the folder of recordings (NeuroLight or Diagnosys, Colordome is for previous machine Astrid excel data) A waiting window appear if you have not yet processed data (3-10s process per recording depending on length frequency and quality)

You see now mean recording with a horizontal slider at index 0 bellow, move the slider (or press -> keyboard) to see individual recordings

You can now validate/refuse artefacts corrections and save data (export/data on menu)

**Group button:** zoom to single flash (a vertical slider appears bind with keyboard -> up/down to change flash)

**Fit button:** to see how the software fits recording (purple curve now, green before)

**Green v button:** valid an artefact detection, select final green curve

**Orange x button\*:** reject an artefact detection, select intermediate orange curve

The orange intermediate curve becomes then final green

**Red X button\*:** reject a flash or whole recording for the one on top left

All treatments disappeared and removed from means/distributions

**Undo all button\*:** go back to initial treatments for the recording you see

**Undo button\*:** undo last click on green v or orange x button.

\* updates plots/means/measures/distributions and take 1-2s (screen freezes)

**Settings** are saved inside the interface, if data have been processed with different settings, you can choose to see previous settings or run and replace with new settings

Settings are advanced and are useful to run faster while developing the method

**Delete threshold:** acceleration threshold to detect and delete points (default: medium)

→ Changes the amount of artefact detected by 1<sup>st</sup> step, misses can be compensated by 2<sup>nd</sup> step but low level may delete reel part of recording

**Drop detection:** minimal size to detect drops and propose artefacts (default: medium)

→ Changes the amount of artefacts detected by 2<sup>nd</sup> step, high level may miss small artefacts and low level may detect too much artefacts and lead to many user rejections

**Fit precision:** maximum time/iteration to fit recordings (default: high)

→ Low level makes the process faster, especially for worst recordings

The precision of the fit is not the most relevant to detect drops

**Derivative smoothing:** smoothing level of Savitzky-Golay filter (default: medium)

→ Uses 8 to 4 polynomial order interpolation for 19 to 27 window points

**See fit:** not plotting fit makes improve a little plot time (default: yes)

**See zoom:** not plotting individual flash makes improve a lot plot time (default: yes)